



Manufactured by **ASTI**

Distributed by:

PromaMedical GmbH
Wernher-von-Braun-Str. 10,
85640 Putzbrunn
Germany
Phone: +49-(0)89-7474-713-0
Fax: +49-(0)89-7474-713-1

PromaMedical Ltd.
The Stables
Egg Lane Hixon
ST18 0PR Staffordshire
United Kingdom
Phone: +49-(0)89-7474-713-0
Fax: +49-(0)89-7474-713-1


PromaMedical Inc
117 N.W. 42nd Avenue, Suite CU-1
Miami, FL, 33126
United States of America
Phone: +1-305-424-1373
Fax: +1-305-460-3529

PromaMedical BV.
Noelweg 22
6171 HB Stein (Maastricht)
Netherlands
Phone: +31-4380-808-49



REVOLUTION OF INJECTION

Length Adjustable
Multi Needle




Length adjustable
Multi-needle attachment



Manufactured by **ASTI**

Special features

- Four injections simultaneously
- High quality needles in 31G or 34G
- Adjustable length of needles up to 4mm
- Insertion at right angle (vertical)



Less injection pressure

The extrusion force was compared between 31G Quatron and 31G conventional single needle connected to 1 ml syringe filled with hyaluronic acid of 12,430 mPa-s.

The plunger was pushed at a rate of 21 μ L/min.



Potential Fields of application

Mesotherapy

Meso
Botulinum Toxin

Local
Anesthesia


Hyaluroic acid
(uncrosslinked)

Lipolysis

and many more!

Advantages

- Less invasive, faster and less downtime
- Less hematoma and less bleeding
- Even and wide distribution of the solution
- Easy change of needle length
- Reduced injection pressure (plunger pressure)



Fewer injections to cover injection site

With Quatron, in comparison with a conventional needle, there is a more even and wider distribution and diffusion of the active substance.